







## Immediate Eyewash Treatment at the Accident Scene AQUA NIT® Sterile Water + ocular NIT® Sterile Solution

Where the eyes or body of any person may be exposed to injurious or corrosive materials/sprays, suitable facilities for quick drenching or flushing of the eyes should be provided within the work area for immediate emergency use.

- Both are medical devices in accordance with European medical device directives
- CE marked in accordance with the German DIN 12 930 temperature withstanding • 3 year shelf life, when unopened • ready for immediate use • no contamination • removal of foreign objects from the eye • can be used single handedly • application by directing the nozzle downwards • have a safety closure controllable liquid emission.

### Eye Wash, AQUA NIT® - Sterile Water

Spray bottle, ready for immediate usage with 250 ml sterile water for an immediate eye wash at the scene of an accident. Suitable for all work environments i.e. in the garage, the factory or working with dangerous materials. Immediate eye wash therapy can prohibit long term serious eye injuries.

20 10 008

## Eye Wash, ocular NIT® - Sterile Solution

Immediate Eye wash sterile solution. 250 ml of sterile solution for eye wash treatment at the accident scene. Especially effective against eye injuries stemming from acidic and caustic burns as well as pepper, CS and tear gas. Three year shelf life when unopened.

20 10 013

#### **AQUA NIT®-Boy Bottle Holder**

Stable metal ring-holder for the secure storage of AQUA NIT® or ocular NIT® eye wash bottles. It can be secured onto a wall or machine by using screws or a double-sided extra strong scotch tape. Convenient and easy to reach storage means fast and effective action at the scene of an accident.

20 10 005

# Eye Wash Stations, AQUA NIT® + ocular NIT®

Eye accidents are some of the most dangerous and common industrial mishaps. These Emergency Eye Wash Stations are the most effective way to deal with eye injuries immediately after the accident. Each box with 4 ready to use spray bottles AQUA NIT® sterile water or ocular NIT® sterile solution, plus first-aid set for eye injuries. Plastic box, size 260 x 160 x 70 mm, orange with printed inscriptions and practical wall brackets. The box can be quickly removed and brought to the scene of an accident scene immediately. The first-aid set for eye injuries contains DuOcul eye compresses, eye flaps, SÖHNGEN® Silk plaster strips, DERMOTEKT® compresses of 10 x 10 cm and disposable gloves.

AQUA NIT® Box with 4 Bottles of AQUA NIT® Sterile water

20 10 009 and Eye Injury Treatment Set

Oculay NIT® Box with 4 Bottles of ocular NIT® Sterile Solution

20 10 004 and Eye Injury Treatment Set

# Eye Wash Bottles, SÖHNGEN® Combi-System

Combination of an eye wash bottle with two replacement bottles. Bath style nozzle attachment for the bathing and cleaning of eyes after light injuries. Tube attachment for pressure cleaning after more serious injuries and the removal of foreign bodies from the eye area. Capacity of 820 ml. With wall bracket.

20 10 006

#### Eye Wash Set, Combi-System

High quality and stylish QUICK-CD case, size 260 x 170 x 110 mm. Manufactured from colour and temperature resistant ABS-plastic material, orange colour, complete with wall bracket and 90° stop-arrest opening system. Two SÖHNGEN®-Combi-Eye-Wash bottles with a capacity of 820 ml.

20 10 014