



Product Specifications

Evolution Plus 2x500ml Eye Wash Kit without Mirror

- Emergency eye irrigation station for any work environment
- Sterile 0.9% saline solution ideal for irrigating eyes
- Meets HSE guidelines for workplaces with a risk of debris or chemical contamination
- Supplied with a quick release wall fixing bracket and kit identification label
- Kit includes: 2x 500ml HypaClens Sterile Eye Wash Bottles and 2x HypaCover Sterile Eye Dressings

The **Evolution Plus 2x 500ml Eye Wash Kit** is a wall-mounted dispenser and ensures that in the event of an emergency, eye wash is easily located. The kit is supplied in a durable hinged translucent plastic Evolution case that includes bottle stabilisers to hold contents firmly in place. The included quick release wall fixing bracket and kit identification label eases storage and identification. Once on the wall bracket, the kit can be opened and worked with as a self-supporting workbench. HypaClens Sterile Eyewash is a single-use medical device intended for emergency eye irrigation. As per HSE guidance, the eye should be washed out with either clean water or sterile fluid from a sealed container for at least 10 minutes following contamination from chemicals or small particles. The eyewash is designed as a substitute for clean running water in first aid situations and is only intended to be used once. Re-use of the bottle would mean you can no longer guarantee the contents are sterile and face the risk of infection. For this reason, the bottle is designed so it cannot be resealed. After use, any excess fluid remaining should be discarded. This eyewash kit assures you have adequate facilities at your place of work. It is available with or without a mirror.



Contents	Qty
HypaClens Eye Wash Bottles, 500ml	2
HypaCover Sterile Eye Dressings	2