



Burn Stop Burns Dressing - 20x20cm

- Gel soaked burn dressings
- Pre cooling not required
- Helps prevent contamination
- Conforms to BS 8599-1 requirements
- Also available in sachets and 50ml bottle size

Burn Stop Burn Dressings provides fast, one-step relief for emergency first aid when someone has been burned. The dressing stops burn progression and cools the skin, relieving pain and protecting against airborne contamination. The gel soaked burn dressing takes no special training to give fast and effective emergency burn care for any type of burn injury, and is available in three sizes depending on your needs. These dressings can also be applied directly to the burn without the need to pre-cool the burn before applying the dressing.

Product Specifications

SKU	D8061
dimensions	20x20cm
dressing_size	20x20cm
brand	burn sstop



TEL: 0845 872 0141